



Rocket® IPC Portfolio is a suite of products and services designed to enable personalized care for people with pleural effusions and malignant ascites.

Rocket® IPC is clinically proven and provides the flexibility to be placed in a procedure suite, an operating room or at the bedside.

- Single catheter - for pleural or peritoneal use
- Complete insertion kit - components for preparing, placing, draining and dressing the catheter
- Drug free option - without lidocaine and skin preparation
- Mini Insertion Kit - VATS option
- Choice of sizes - dilators, echogenic needles

Code	Product	Catheter	Kit Contents		
R51400-MI-PL	Rocket® IPC™ Mini Insertion Kit (VATS)	<i>Rocket® IPC™ Silicone Catheter 16FG with metal tunneller</i>	<ul style="list-style-type: none"> • 1 x drainage line • 1 x valve cap • 1 x tunnelling forceps • 1 x foam catheter pad • 1 x clear-film dressing 	MINI KIT	
R51400-16-00 <i>Includes all items contained in R51400-MI-PL</i>	Rocket® IPC™ Insertion Kit - Drug Free	<i>Rocket® IPC™ Silicone Catheter 16FG with metal tunneller</i>	<ul style="list-style-type: none"> • 1 x 18G Echogenic Guidewire insertion needle (70mm) • 1 x 18G Echogenic Guidewire insertion needle (100mm) • 1 x J-Tip Guidewire with markings (70 cm) • 1 x 16 Fr Split Sheath / Dilator (11 cm) • 1 x 16 Fr Split Sheath / Dilator (17 cm) • 1 x 14 Fr Dilator 8" • 1 x 14 Fr Dilator 10 cm • 1 x IPC Catheter Stiffener • 1 x No. 11 Blade Safety Scalpel • 1 x 10ml syringe • 1 x 21G x 1.5" Green Safety Needle • 1 x 25G 5/8" Orange Safety Needle • 1 x Absorbable Suture • 1 x 2.0 Suture Curved Needle • 1 x Needle safety vise • 1 x 2 part Adhesive Fenestrated Drape • 1 x 2 part Large patient drape • 1 x Single use Suture Forcep • 1 x 18G Blunt fill needle and filter • 1 x IPC insertion label for chart • 1 x pair Latex-free Gloves • 8 x Gauze pads • 1 x Scissors • 1 each - Wide bore adapter, Extension tubing, Needleless port adapter • 1 x Tweezers 	DRUG FREE KIT	COMPLETE INSERTION KIT
R51400-16- PL <i>Includes all items contained in R51400-16-00</i>	Rocket® IPC™ Complete Insertion Kit	<i>Rocket® IPC™ Silicone Catheter 16FG with metal tunneller</i>	<ul style="list-style-type: none"> • Lidocaine (1%, 5 mL x 3) • Skin preparation (6 mL x 2) 		

Code	Product Description	Code	Product Description
R54410-00-CP	Rocket® IPC™ Valve Cap Sterile, for single use, 10 per box	R54410-WB-DL	Rocket® IPC™ Wide Bore Drainage Line Sterile for single use, 10 per box
R54410-00-VV	Rocket® IPC™ Replacement Valve Kit Supplied individually, sterile for single use, 10 per box	R54410-00-DL	Rocket® IPC™ Drainage Line Sterile for single use, 10 per box
R54400-00-RK	Rocket® IPC™ Catheter Removal Kit Supplied individually, sterile for single use, 5 per box	R54410-16-11	Rocket® IPC™ Split Sheath 16G x 11cm Sterile for single use, 10 per box
R54400-00-00	Rocket® IPC™ Dressing Pack Sterile for single use, 5 per box	R54410-16-17	Rocket® IPC™ Split Sheath 16G x 17cm Sterile for single use, 10 per box

