

Rocket® Seldinger Chest Drainage Sets



The **Rocket® Seldinger Chest Drainage Set** is indicated for the percutaneous insertion of a pleural drain for fluid and air and is particularly suitable for: ¹

- Pneumothorax, spontaneous or iatrogenic
- Effusions, benign or malignant, except those with high protein content.
- Low viscosity empyema

The **Rocket® Seldinger Chest Drainage Set** features:

- Improved patient comfort using small bore catheterisation ^{2,3}
- Reliable insertion via Seldinger approach under direct vision, x-ray screening or ultrasound guidance
- 10cm dilators with safety stops to minimise the risks of over-insertion associated with the traditional 20cm dilator.
- Drapes and a gown are included in all PK coded packs in response to NPSA/2008/RRR3

References:

1. *Clinical Evaluation of Seldinger Drainage Kit*. Beaumont A C et al Proceedings of Intensive Care Society Meeting. December 1999.
2. *Catheters smaller than 24 French gauge can be used for chest drains*. Taylor PM. BMJ (1997) 315:186
3. *Management of Malignant Pleural Effusions* Antunes G, Neville

CODE	CODE	CODE
Rocket® Seldinger Chest Drainage Set with Safety Guards Skin knife, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir-tree & wide-bore luer connectors. Sterile, single use in cartons of 5 units.	12FG	R54544-12-SG USL 854544-12-SG
	18FG	R54544-18-SG USL 854544-18-SG
Rocket® Seldinger Chest Drainage Tray Complete Ward Procedure Set Skin knife, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir-tree & wide-bore luer connectors, skin fixation set, ward procedure pack including drapes, swabs, sponges, 10ml syringes, gallipot. Sterile, for single use in cartons of 5 units.	12FG	R54544-12-PK USL 854544-12-PK
	18FG	R54544-18-PK USL 854544-18-PK
	20Fg	R54549-20-PK USL 854549-12-PK